



Material-friendly rapid cleaning and disinfection wipes impregnated with an low-alcoholic solution.

mikrozid® universal wipes premium

Our Plus

- good material compatibility (e.g. tablets)
- rapid effect (Noro 30 sec)
- effective according to EN 16615 (4-field test) in 1 minute
- dermatologically tested
- high quality and big wipes
- colourant and fragrance free
- outstanding cleaning performance

Application areas

Alcohol-based rapid disinfection of medical devices in all areas and all kinds of wipeable surfaces with an increased risk of infection, where short contact times are required and where a material friendly disinfectant is needed, e.g.:

- Patient nearby surfaces
- Patient treatment units
- Examination couches
- Operating tables and nearby work surfaces
- Surfaces of medical equipment
- Keyboards or touch pads of sensitive communication equipment, e.g. smartphones or tablets

Instructions for use

Thoroughly wipe the surface with the impregnated cloth and allow to take effect. Make sure to wet surfaces completely and keep them wet for the whole exposure time. Make sure that all visible dirt is removed prior to disinfection. When opening the softpack, make sure that hygiene requirements are maintained. Make sure to close the lid tightly after use. Not for final disinfection of semicritical and critical medical devices! schülke guarantees a shelf life of 1 month after the package has been opened. The product can be used without gloves, if allowed by infection and occupational protection.

Microbiological efficacy

Efficacy	Concentration	Contact time
bactericidal EN13727 - dirty conditions	ready-to-use	15 sec.
tuberculocidal EN14348 - dirty conditions	ready-to-use	5 min.
yeastcidal EN13624 - dirty conditions	ready-to-use	1 min.
virucidal against enveloped viruses - dirty conditions	ready-to-use	15 sec.
limited spectrum virucidal activity EN14476 - dirty conditions	ready-to-use	15 min.
Norovirus EN14476 - clean conditions	ready-to-use	30 sec.
Norovirus EN14476 - dirty conditions	ready-to-use	60 sec.
Polyoma SV40 - dirty conditions	ready-to-use	15 sec.
Rotavirus EN14476 - clean conditions	ready-to-use	15 sec.

Certificates

- IHO listed
- VAH certificate



Product data

Composition:

One softpack contains 100 ready-to-use pre-saturated wipes.
100 g solution contains the following active substances: 17,4 g Propan-2-ol, 12,6 g Ethanol (94 % w/w)
Labelling according to Regulation (EC) No. 648/2004: < 5% anionic surfactants

Chemical-physical data

Color	colourless
Density	ca. 0,95 g/cm ³ / of the active solution
Flash point	26 °C / Method : DIN 51755 Part 1 / of the active solution
Form	Aqueous, alcohol containing solution on non-woven
pH	3 - 3,6 / of the active solution
Viscosity, dynamic	No data available

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining out high quality standards.

Expert opinion and information

Please visit our website for an overview of all available literature/reports on the product: www.schuelke.com.

For individual questions:

Customer Sales Service

Phone: +49 40 52100-666

E-Mail: info@schuelke.com

*The surface coverage of one wipe depends on the surrounding conditions and the structure of the wiped surface.

Special advice

Use disinfectants safely. Always read the label and product information before use.

The product is highly compatible with metals and synthetics. Please follow carefully any national fire and explosion protection guidelines when using alcohol-based disinfectants. Only treat small surfaces. The wiping range depends on the room temperature and on the structure of the surface to be disinfected. More information available on request.

Following sizes are available:

- mikrozid® universal wipes *premium* (6 x 100 wipes 20 x 20 cm) Surface coverage*: approx. 1.0 m²
- mikrozid® universal wipes *premium maxi* (6 x 80 wipes 25 x 25 cm) Surface coverage*: approx. 1.5 m²

Information for order

Item	Delivery form	Item no.
mikrozid® universal wipes premium (100 wipes)	6/Carton	on request
mikrozid® universal wipes premium maxi (80 wipes)	6/Carton	on request

These products are not available in every country. For more information please contact our local subsidiary or distributor.

Related Products

- mikrozid® AF wipes
- mikrozid® PAA wipes
- mikrozid® sensitive wipes
- mikrozid® universal liquid



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

Schülke Headquarters
Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 40 - 52100 - 0
Fax +49 40 - 52100 - 318
www.schuelke.com
info@schuelke.com