

**SARS-CoV-2 Antigen Rapid Test Kit**  
**(Fluorescence Immunochromatography)**

**Clinical Evaluation Report**

**30 November 2020**

*Collaboration:*

*Wiesbaden Germany, Microbiology Institute of German Armed Forces*

## Objective

The objective of this study was to evaluate the diagnostic performances of SARS-CoV-2 Antigen Rapid Test Kit (fluorescence immunochromatography) produced by Biohit Healthcare (Hefei) Co., Ltd. (sensitivity, specificity and positive and negative predictive values from pre-specified prevalence assumptions) carried out on nasopharyngeal samples collected.

## Material and methods

This was a retrospective study carried out using nasopharyngeal samples, from 119 patients with COVID-19 infection confirmed by PCR. Negative samples were collected from 250 patients with confirmed PCR negative results.

Samples tested included:

- **119** nasopharyngeal samples frozen at -80 ° C in 3ml NaCl or 3ml PBS, **tested positive for SARS-CoV-2 by PCR** (RT-PCR by Roche Cobas system), tested in November, 2020. The time span of the samples is from beginning of the pandemic to November 2020. The **70** samples from **Labor Dr. Riegel in Wiesbaden Germany**, and **49** samples from **Microbiology Institute of German Armed Forces**. The nasopharyngeal swab samples were collected on the day of confirmed diagnosis of COVID-19 infection (the day of PCR test result is positive).
- **250** nasopharyngeal samples frozen at -80 ° C in 3ml NaCl or 3ml PBS from the same sources, **tested negative for SARS-CoV-2 by PCR** (RT-PCR by Roche Cobas system), tested in November.

Reagent and Instrument:

- SARS-CoV-2 Antigen Rapid Test Kit (fluorescence immunochromatography)  
Manufacturer: Biohit Healthcare (Hefei) Co., Ltd.; Batch number: SA200903.
- Ultraviolet flashlight: TANK007 (365nm wavelength)  
Manufacturer: Shenzhen Zhong Dao Electronics Co., Ltd.

## Results

The Biohit results are shown below table:

- 1) **Sensitivity** (Se, percentage of positive RAPID results among the cases identified as positive by PCR [i.e., true positives / (true positives + false negatives)]);
- 2) **Specificity** (Sp, percentage of negative RAPID results among the cases identified as negative [i.e., true negatives / (true negatives + false positives)]);
- 3) **Positive predictive values** (PPV: probability that a positive test is a true positive in PCR, defined by the formula  $Se * P / (Se * P + (1-P) * (1-Sp))$  [with P = prevalence]); calculated for pre-specified prevalence of infection in the tested population of 1% and 5%.

4) **Negative predictive values** (NPV: probability that a negative test is a true negative in PCR, defined by the formula  $Sp * (1-P) / (Sp * (1-P) + P * (1-Se))$  [with P = prevalence]); calculated for pre-specified prevalence of infection in the tested population of 1% and 5%.

Specificity				Sensitivity				PPV		NPN	
N	%Negative	IC 95%		N	%Positive	IC 95%		P=1%	P=5%	P=1%	P=5%
250	99.6%	97.79%	99.99%	119	78.15%	69.65%	85.20%	66.37%	91.14%	99.78%	99.77%
CT Value				Number	Sensitivity	IC 95%					
CT ≤ 20				8	100%	63.06%	100.00%	71.63%	92.94%	100.00%	100.00%
CT ]20-25]				36	100%	90.26%	100.00%	71.63%	92.94%	100.00%	100.00%
CT ]25-30]				46	84.78%	71.13%	93.66%	68.16%	91.77%	99.85%	99.84%
CT > 30				29	34.48%	17.94%	54.33%	46.54%	81.94%	99.34%	99.31%
CT ≤ 33				106	86.79%	78.83%	92.59%	68.67%	91.95%	99.87%	99.86%
CT ≤ 30				90	92.22%	84.63%	96.82%	69.96%	92.39%	99.92%	99.92%
CT ≤ 25				44	100%	91.96%	100.00%	71.63%	92.94%	100.00%	100.00%
<b>Overall sensitivity versus PCR: 78.15% - Sensitivity for Ct ≤ 33: 86.79% - Specificity: 99.6%</b>											

### Comparison with other COVID- 19 antigenic TRODs

Source : Diagnostic performance evaluation of covid-19 rapid diagnostic localization test (<https://www.aphp.fr/contenu/evaluation-de-la-performance-diagnostique-de-s-tests-rapides-dorientation-diagnostique>)

This was a retrospective study carried out using nasopharyngeal samples collected prospectively. Samples tested included **297 aliquots** of nasopharyngeal samples frozen at -80 ° C in commercial viral transport medium (Cepheid® or Deltalab®) or physiological serum, **tested positive for SARS-CoV-2 by PCR, and 337 aliquots** of nasopharyngeal samples frozen at -80 ° C in commercial viral transport medium (Cepheid®) **negative for SARS-CoV-2** because they were taken before the virus circulating period, i.e., between April and August 2019. Six COVID-19 antigenic TRODs were verified, **the ABBOTT, BIOSYNEX and AAZ tests were the most efficient.**

## 1) Test ABBOTT

ABBOTT												
	Spécificité				Sensibilité				V P		V P N	
	N témoins	Témoins négatifs	IC 95%		N cas	% C as positifs	IC 95%		Prévalence 1%	Prévalence 5%	Prévalence 1%	Prévalence 5%
<b>Global</b>	337	100,0%	98,9%	100,0%	295	<b>55,3%</b>	49,4%	61,0%	100,0%	100,0%	99,6%	97,7%
<b>Par sous-groupes</b>	Délai d'apparition des symptômes											
		Délai 0-3j			97	<b>79,4%</b>	70,0%	86,9%	100,0%	100,0%	99,8%	98,9%
		Délai 4-7j			103	<b>52,4%</b>	42,4%	62,4%	100,0%	100,0%	99,5%	97,6%
		Délai 8-11j			63	<b>33,3%</b>	22,0%	46,3%	100,0%	100,0%	99,3%	96,6%
		Délai ≥12j			24	<b>37,5%</b>	18,8%	59,4%	100,0%	100,0%	99,4%	96,8%
		Délai ≤7j			200	<b>65,5%</b>	58,5%	72,1%	100,0%	100,0%	99,7%	98,2%
		Ct value										
		Ct ≤20			40	<b>95,0%</b>	83,1%	99,4%	100,0%	100,0%	99,9%	99,7%
		Ct ]20-25]			90	<b>83,3%</b>	74,0%	90,4%	100,0%	100,0%	99,8%	99,1%
		Ct ]25-30]			73	<b>57,5%</b>	45,4%	69,0%	100,0%	100,0%	99,6%	97,8%
		Ct >30			88	<b>8,0%</b>	3,3%	15,7%	100,0%	100,0%	99,1%	95,4%
		Ct ≤33			245	<b>65,7%</b>	59,4%	71,6%	100,0%	100,0%	99,7%	98,2%
		Ct ≤25			130	<b>86,9%</b>	79,9%	92,2%	100,0%	100,0%	99,9%	99,3%
		Ct ≤23			96	<b>94,8%</b>	88,3%	98,3%	100,0%	100,0%	99,9%	99,7%
	Sévérité											
	Bénin			202	<b>58,4%</b>	51,3%	65,3%	100,0%	100,0%	99,6%	97,9%	
	Sévère			92	<b>47,8%</b>	37,3%	58,5%	100,0%	100,0%	99,5%	97,3%	

Sensibilité globale *versus* PCR : **55,3%** - Sensibilité pour Ct ≤33 : **65,7%** - Spécificité : **100%**  
 NB : Résultats invalides : 2 positifs en PCR

Overall sensitivity *versus* PCR: 55.3% - Sensitivity for Ct ≤33: 65.7% - Specificity: 100%  
 NB: Invalid results: 2 positive in PCR

## 2) BIOSYNEX

BIOSYNEX												
	Spécificité				Sensibilité				V P		V P N	
	N témoins	Témoins négatifs	IC 95%		N cas	% C as positifs	IC 95%		Prévalence 1%	Prévalence 5%	Prévalence 1%	Prévalence 5%
<b>Global</b>	337	98,5%	96,6%	99,5%	297	<b>59,6%</b>	53,8%	65,2%	28,6%	67,6%	99,6%	97,9%
<b>Par sous-groupes</b>	Délai d'apparition des symptômes											
		Délai 0-3j			97	<b>81,4%</b>	72,3%	88,6%	35,4%	74,1%	99,8%	99,0%
		Délai 4-7j			103	<b>56,3%</b>	46,2%	66,1%	27,5%	66,4%	99,6%	97,7%
		Délai 8-11j			63	<b>42,9%</b>	30,5%	56,0%	22,4%	60,1%	99,4%	97,0%
		Délai ≥12j			26	<b>42,3%</b>	23,4%	63,1%	22,2%	59,8%	99,4%	97,0%
		Délai ≤7j			200	<b>68,5%</b>	61,6%	74,9%	31,6%	70,6%	99,7%	98,3%
		Ct value										
		Ct ≤20			40	<b>97,5%</b>	86,8%	99,9%	39,6%	77,4%	100,0%	99,9%
		Ct ]20-25]			90	<b>92,2%</b>	84,6%	96,8%	38,3%	76,4%	99,9%	99,6%
		Ct ]25-30]			74	<b>63,5%</b>	51,5%	74,4%	30,0%	69,0%	99,6%	98,1%
		Ct >30			89	<b>9,0%</b>	4,0%	16,9%	5,7%	24,0%	99,1%	95,4%
		Ct ≤33			247	<b>71,3%</b>	65,2%	76,8%	32,4%	71,4%	99,7%	98,5%
		Ct ≤25			130	<b>93,8%</b>	88,2%	97,3%	38,7%	76,7%	99,9%	99,7%
		Ct ≤23			96	<b>96,9%</b>	91,1%	99,4%	39,5%	77,3%	100,0%	99,8%
	Sévérité											
	Bénin			202	<b>60,4%</b>	53,3%	67,2%	28,9%	67,9%	99,6%	97,9%	
	Sévère			94	<b>58,5%</b>	47,9%	68,6%	28,3%	67,2%	99,6%	97,8%	

Sensibilité globale *versus* PCR : **59,6%** - Sensibilité pour Ct ≤33 : **71,3%** - Spécificité : **98,5%**  
 NB : Résultats invalides : aucun

Overall sensitivity *versus* PCR: 59.6% - Sensitivity for Ct ≤33: 71.3% - Specificity: 98.5%  
 NB: Invalid results: none

### 3) AAZ

AAZ												
	Spécificité			Sensibilité				V P		V P N		
	N témoins	Témoins négatifs	IC95%	N cas	% C as positifs	IC95%		Prévalence 1%	Prévalence 5%	Prévalence 1%	Prévalence 5%	
<b>Global</b>	337	100,0%	98,9%	100,0%	295	<b>61,7%</b>	55,9%	67,3%	100,0%	100,0%	99,6%	98,0%
Délai d'apparition des symptômes												
Délai 0-3j				97	<b>81,4%</b>	72,3%	88,6%	100,0%	100,0%	99,8%	99,0%	
Délai 4-7j				103	<b>61,2%</b>	51,1%	70,6%	100,0%	100,0%	99,6%	98,0%	
Délai 8-11j				63	<b>42,9%</b>	30,5%	56,0%	100,0%	100,0%	99,4%	97,1%	
Délai ≥12j				24	<b>37,5%</b>	18,8%	59,4%	100,0%	100,0%	99,4%	96,8%	
Délai ≤7j				200	<b>71,0%</b>	64,2%	77,2%	100,0%	100,0%	99,7%	98,5%	
Ct value												
Ct ≤20				40	<b>100,0%</b>	91,2%	100,0%	100,0%	100,0%	100,0%	100,0%	
Ct ]20-25]				90	<b>94,4%</b>	87,5%	98,2%	100,0%	100,0%	99,9%	99,7%	
Ct ]25-30]				73	<b>65,8%</b>	53,7%	76,5%	100,0%	100,0%	99,7%	98,2%	
Ct >30				88	<b>9,1%</b>	4,0%	17,1%	100,0%	100,0%	99,1%	95,4%	
Ct ≤33				245	<b>73,5%</b>	67,5%	78,9%	100,0%	100,0%	99,7%	98,6%	
Ct ≤25				130	<b>96,2%</b>	91,3%	98,7%	100,0%	100,0%	100,0%	99,8%	
Ct ≤23				96	<b>97,9%</b>	92,7%	99,7%	100,0%	100,0%	100,0%	99,9%	
Sévérité												
Bénin				202	<b>62,4%</b>	55,3%	69,1%	100,0%	100,0%	99,6%	98,1%	
Sévère				92	<b>59,8%</b>	49,0%	69,9%	100,0%	100,0%	99,6%	97,9%	

Sensibilité globale versus PCR : **61,7%** - Sensibilité pour Ct ≤33 : **73,5%** - Spécificité : **100%**  
 NB : Résultats invalides : 2 positifs en PCR

Overall sensitivity versus PCR: **61.7%** - Sensitivity for Ct ≤33: **73.5%** - Specificity: **100%**  
 NB: Invalid results: 2 positive in PCR

### Conclusion

Based on the above research, Biohit SARS-CoV-2 Antigen Rapid Test Kit (fluorescence immunochromatography) is most effective in diagnostic performances of COVID-19.