



CE
certification

COVID-19 Antigen Rapid Test Cassette (Nasal Swab)






- Designed for qualitative detection of SARS-CoV-2 antigen in samples

■ Antigen rapid tests


Ag rapid tests are useful for detecting SARS-CoV-2 antigens, most often nucleocapsid, produced by the replicating virus in respiratory secretions.

Based on the above principle, COVID-19 Antigen Rapid Test Cassette (Nasal Swab) manufactured by Clongene is intended for the qualitative detection SARS-CoV-2 nucleocapsid antigens from individuals who are suspected of COVID-19 by their healthcare provider.

■ Kit Contents


Work Station	Extraction Reagent	Test Cassette	Extraction Tube & Dropper Tip	Sterilized Swab
				

■ Principle




■ Procedure

sample collection




Insert swab about 2.5 cm (1 inch) into nostril until resistance is met at turbinates and rotate the swab.

test procedure



■ Features



Hangzhou Clongene Biotech Co., Ltd.

<http://en.clongene.com/>

E-mail: marketing@clongene.com

Contact: +86-571-88262120

