

Instruction for Use Urine bag		Prepared by	Li Xiuren
		Checked by	Zhou Ting
Doc. No.	KB/CE01-08	Approved by	Zhou Hongbao
Effective date	2022.03.01	Ver. B/0	Page No. Page 1 of 2

Instruction for Use

[Product name]: Urine bag

[Specification]: 1000ml, 1500ml, 2000ml, 2000+200ml, 2000+500ml, 2600+400ml, 2600+500ml.

[Materials]: PVC film (bag), PVC (connecting tube, connector), PP (hook etc.)

[General Device Description]: The urine bag is composed of Bag, cross valve seat, Cross valve, liquid inlet seat, breathe film, Clappet, Hook, Bandage, liquid inlet tube, Sheets clamp, sample connector, protective cap.

[Intended use]: The urine bag connecting with drainage tube, used in storing the urine, For Clinical Use. Non direct contact with the body. For single-use. Use time: short term.

[Sterilization Method]: EO sterilization

[Validity]: 5 years.

[Use method]:

- Check the connection between the bag and the connecting tube.
- Connect the urine bag with nelaton tube, settle down the urine bag.

[Cautions]:

- For single patient use only.
- Sterile if package is unopened or damaged.
- Store in a dry, cool and dark place.

[Storage and transportation]:

1. The device shall be stored at a relative humidity of less than 80%, the storage environment shall be kept non-corrosive gases, dry, avoiding sunlight, well-ventilated conditions, and clean.

2. The device shall be taken care during transportation and handling.

[Contact information]:

Changshu Kangbao Medical Appliance Factory

Furong Village, Baimao, Guli Town, 215500 Changshu City, Jiangsu Province, PEOPLE'S REPUBLIC OF CHINA

TEL: +86-512-52539293 FAX: +86-512-52301322

Email: kb@kang-bao.cn Website: www.kang-bao.cn

The authorized representative of European Union:

MedPath GmbH

Mies-van-der-Rohe-Strasse 8,

80807 Munich, Germany

Phone: +49(0)89 189174474

Fax: +49(0)89 5485 8884

Email: info@medpath.pro

[Symbols]

Instruction for Use Urine bag			Prepared by	Li Xiuren
			Checked by	Zhou Ting
Doc. No.	KB/CE01-08		Approved by	Zhou Hongbao
Effective date	2022.03.01	Ver. B/0	Page No.	Page 2 of 2

	CE Mark		Lot number
	Single use		Date of manufacture
	Sterilized by EO		Manufacturer
	Notice		European union representative
	The product includes DEHP		Indicates the item is a medical device